# Gesamtkatalog Labor

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Carbide cutters for dental CAD/CAM systems

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Recommendations for use and safe operation

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stands for high performance carbide instruments from acurata, Germany. Manufactured from micro-grain tungsten carbide in state-of-art CNC machines, these instruments meet the highest demands of technicians. Precisely milled blade geometries, high efficiency and vibration-free concentricity ensure sparing drives and materials. Quality is managed through all manufacturing stages according to DIN EN ISO 13485. These instruments are characterised by high performance and efficiency in the daily laboratory work.

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

The toothings of **acurata** cutters for laboratory fulfil the different requirements for optimal grinding surface patterns.

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


## Toothng-navigator . Here you find the appropriate cutter for your material

Material / work step	toothng	recommended speeds	colour code
<b>Moulding plaster</b>			
Coarse material removal on wet plaster	223 / SGF	10.000 min <sup>-1</sup>	
Coarse material removal on dry plaster	220	12.000 min <sup>-1</sup>	
Work on model stumps	175 / 190	12.000 min <sup>-1</sup>	
<b>Precious metals</b>			
Inlays, Onlays, Bridges	190 / RRX	15–20.000 min <sup>-1</sup>	
Elaborate, recontour	140 / RRX	15–20.000 min <sup>-1</sup>	
Adjustment, abrade surfaces	140 / RRX	15–20.000 min <sup>-1</sup>	
Elaborate chewing surfaces	110	10–15.000 min <sup>-1</sup>	
Roughen of surfaces which are to be blinded			
<b>Model casting</b>			
Basis of prosthesis	190 / 134 / NEM / RRX	15.000 min <sup>-1</sup>	
Elaborate, recontour	140 / NEM / RRX	15.000 min <sup>-1</sup>	
Adjustment, abrade surfaces			
<b>None precious alloys</b>			
Inlays, bridges	190 / 134 / NEM / RRX	15.000 min <sup>-1</sup>	
Elaborate recontour	140 / 134 / NEM / RRX	15.000 min <sup>-1</sup>	
Adjustment, abrade surfaces	140 / NEM / RRX	15.000 min <sup>-1</sup>	
Elaborate chewing surfaces	110 / NEM / RRX	12.000 min <sup>-1</sup>	
Roughen of surfaces which are to be blinded			
<b>Titan, cast</b>			
Crowns, bridges, basis of prosthesis	154	15.000 min <sup>-1</sup>	
Elaborate, recontour			
<b>Fibre-reinforced acrylics</b>			
Moulding material	176	12.000 min <sup>-1</sup>	
Elaborate			
Crowns, Bridges	134	15.000 min <sup>-1</sup>	
Elaborate, recontour			
<b>Facing acrylics</b>			
Facings	176 / 190 / 134	15.000 min <sup>-1</sup>	
Elaborate, recontour	140	15.000 min <sup>-1</sup>	
Adjustments, abrade surfaces			
<b>Temporary materials</b>			
Temporary appliance	190 / 134	15.000 min <sup>-1</sup>	
Elaborate, abrade			
<b>Denture resin</b>			
Full-/ partial denture	176 / 175 / 190	15.000 min <sup>-1</sup>	
Elaborate			
Soft acrylics	176	15.000 min <sup>-1</sup>	
Elaborate	LGQ	15.000 min <sup>-1</sup>	
<b>Ceramics</b>			
Elaborate chewing surfaces	110 / 133	15.000 min <sup>-1</sup>	
<b>Titan, milled</b>			
	110 / 190 / RRX	15.000 min <sup>-1</sup>	

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### Tungsten carbide cutter with ac-blue wear protection coating

- cutting edge protection for longer service life
- high heat resistance e.g. during dry machining
- universal applicability for all millable materials
- very thin coating thickness of  $1.5 \mu\text{m}$  provides the usual sharp cutting edges
- reduced heat development for a sensitive processing
- simplified cleaning

ac blue provides best results in many different application areas and at the same time maximum efficiency.

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High-quality milling tools for the **treatment of all standard materials** used in the dental industry – from our own production “made in Germany”.

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## RRX-toothing

is the further development of the successful NEM-tothing. Universal cutter for all metals, resin and composite

- coarse and fine segments on one instrument
  - german Patent registered – DE 10 2007 009 304.9
  - **Applied for EU Patent - EP 1 958 589 B1**

The special staggered toothing with a right-wound transversal cut creates very short and sharp cutting edges. Short chips are formed during the milling. These small chips are removed easily and grinded away during the rotation of the cutter. Even with hard materials no needle-shaped chips are generated which could unpleasantly irritate the skin. With the RRX-toothing you achieve smooth and easy polishable surfaces on all metals, resin and composite. The efficient cutting performance, operating comfort and long service life will convince you.

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**Quickly and effectively** – quick removal, smooth running and long service life make our NEM-instruments suitable for **universal applications**.

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TC-Cutter for pre- and superfine milling of non-precious alloys.

The instruments have 2 different cuts (coarse and fine on one tool!) and two chip-breaker grooves for fine- and coarse-grained abrasion. The sharp blades penetrate easily and quickly even into the hardest none-precious metals and CrCo alloys. Due to the different cut structure, the “combination cutters” produce good material removal and smooth, fine surfaces.

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Deutsches Patent angemeldet – DE 10 2007 009 304.9

**EU Patent – EP 1 958 589 B1**

2 different cuts on one tool. Applied for German Patent – DE 10 2007 009 304.9

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**TC ball nose cutter for all alloys.**

The instruments have a new special toothing with cross-cut for short, grainy chips. The sharp cutting edges require only a low contact pressure. They penetrate easily and quickly even the hardest alloys. These instruments complete our NEM-toothing (patent De 10 2007 009 304 9) which is proved over the years. Easily and clearly to identify by their distinctive, orange colour ring.

**Spezialverzahnung mit Querhieb – 8 Schneiden bei allen Durchmessern**  
Special toothing with cross cut – all diameters have 8 blades  
denture spéciale hélicoïdale – 8 lames chez tous les diamètres  
Dentatura con taglio trasversale – 8 tagli presso tutti i diametri

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**Recommendations:** Optimum Speed: 20.000–25.000 min<sup>-1</sup>

with light contact pressure. Max. Speed: 50.000 min<sup>-1</sup>

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Twisted cut, supercoarse with guide phase and s-shaped central cutting edge for plaster

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Lamellas toothing for preparing of prostheses-, interim prostheses and deep-drawing film acrylics

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Twisted cut, fine for fissures, all precious alloys and ceramics

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Twisted cut, medium for plaster and acrylics

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Twisted cut with fine cross cut for gold-,  
Pd.-alloy, none-precious alloys and acrylics

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Twisted cut with cross cut for acrylics

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For titanium casts.

acurata has simplified the handling of titanium casts. Specially designed cutting edge geometry and cut phases with fine cross cut in the acurata cutter for titanium join with the metallurgic attributes of the titanium of the different gradients. The casts remain grooveless, the density of the material unaffected, the crystal lattice for the laser sealing is preserved. The attributes which are, as is generally known, very biocompatible remain untouched.

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Double cut, fine – for all alloys and solid acrylics

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Double cut, medium – for Cr Co, precious alloys and solid acrylics

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Our high performance carbide instruments, specially developed for the dental laboratory, assist dental technicians with their daily work. Our professional, time-saving first class tools for economical work in the dental laboratory.

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Double cut, coarse for plaster and acrylics

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Double cut, supercoarse for plaster and acrylics

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### TC-cutter for fissures

Twist cut to the right and cross cut for standard applications.

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

Twisted cut, to the right with fine cross cut

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### RQF cutter

fine spiral toothing with fine cross cut

- Pre- and fine milling of non-precious alloys and gold alloys with one cutter
- by means of speed and contact pressure, the intended characteristic can be realized
- depending on the material / alloy mostly in the range 20.000–25.000 min<sup>-1</sup>
- fast cutting and smooth running
- preferable to use milling model

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**XNE / FNE - toothing**

for hardly machinable non precious alloys, measurable larger removal performance, smooth, plain surface, unequalled endurance

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**roughing cutter**

high cutting performance, e.g. abutment

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**Parallel TC-Cutter**

round, handpiece, short Ø 2,35 mm, left-hand twist, right cutting.

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**Taper TC-Cutter**

round, handpiece, short Ø 2,35 mm, left-hand twist, right cutting

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**Parallel- and taper- TC-Cutter**

for pure titanium cast. Pre-cutter, left-hand twist, right cutting, shank Ø 2,35 mm.

Micro-milling for titanium we recommend cut 103 (white ring)!

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


### Carbide cutters for dental CAD/CAM systems

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On the following pages, you will find a selection of our most popular CAD/CAM tools. You can find the complete product range either online [www.acurata.de](http://www.acurata.de) or in our CAD/CAM Special folder. Do you have questions?  
Simply call us under phone number +49(0)8504-9117-68.

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Tungsten carbide, ac-blue coating, 6 mm shank diameter, 48 mm allover length, cylindric radius cutter, 2 blades, length 15,0 mm, diameter of the working part 3,0 mm

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For tools of the universal range. The letter code A after the ninth digit = on stock with stop ring. Upon customer request, all tools can be provided with "Vario"-stop ring for 6 mm shank diameter or plastic ring for 3 mm shank diameter.

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Tool geometry as well as clamping length compatible with Amann/Girrbach-cutter

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Tool geometry as well as clamping length compatible with ZirkonZahn-cutter

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Tool geometry as well as clamping length compatible with VHF-cutter

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**1-blade, round**


exclusively for PMMA, PEEK and wax, **shank Ø 3 mm**

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**2-blades, round**

exclusively for ZrO, Sintron, PMMA and PEEK, **shank Ø 3 mm**

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
2-blades, round, metal treatment  
exclusively for CrCo, with spacer ring,  
**shank Ø 3 mm**

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
**2-blades, round**


exclusively for ZrO, PMMA and PEEK,  
with spacer ring, **shank Ø 3 mm**

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**2-blades, round**

**ac-crystalline coated, VHF**

exclusively for ZrO, with spacer ring, **shank Ø 3 mm**

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**1-blade, round**

exclusively for PMMA, PEEK and wax, with spacer ring, **shank Ø 3 mm**

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


"Lollipop" ball cutter,  
with spacer ring  
Milling of undercuts,  
for all milling machines, **shank 3 mm**

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

2- and 3-bladed tools, round, ac-blue coating  
exclusively for ZrO, PMMA and PEEK, with spacer ring, **shank Ø 3 mm**

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



**2-blades, round**


exclusively for ZrO, PMMA and PEEK, with spacer ring  
compatible with "Wieland Select", VHF cam 5 S1 + S2, **shank Ø 3 mm**

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

**2 blades round**


ac-crystalline coated  
exclusively for ZrO, with spacer ring, **shank Ø 3 mm**

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
**2-blades, round, ac-blue-coating**


exclusively for ZrO, PMMA and PEEK, **shank Ø 3 mm**

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2-blades, round, ac-crystalline-diamond coating  
 exclusively for ZrO, **shank Ø 3 mm**

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
2-blades, round  
 exclusively for ZrO and PMMA, **shank Ø 6 mm**

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
2-blades, flat, face cutter  
exclusively for CrCo, **shank Ø 6 mm**

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
2-blades, round  
exclusively for CrCo, **shank Ø 6 mm**

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
2-blades, round,  
with spacer ring, long version  
exclusively for ZrO and PMMA,  
**shank Ø 3 mm until blank height 25 mm**

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
2-blades, conical, Treatment of fissures,  
with spacer ring, long version  
exclusively for ZrO and PMMA, **shank Ø 3 mm until blank height 25 mm**

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
2-blades, round, ac-crystalline diamond coating,  
with spacer ring, long version  
exclusively for ZrO, **shank Ø 3 mm until blank height 25 mm**

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2-blades, conical, round  
exclusively for ZrO, PMMA and PEEK, **shank Ø 6 mm**,  
also available with/without spacer ring

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### 2-blades, round

exclusively for ZrO, PMMA and PEEK, **shank Ø 6 mm**,  
also available with/without spacer ring

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### 2-blades, round

ac-crystalline diamond coating  
exclusively for ZrO, **shank Ø 6 mm**, also available with/without spacer ring

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


Compatible with I-MES, Roeders, Wissner, MB, Primacon, DMG, Dental Concept Systems a.o.

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
2-blades, round, ac-crystalline diamond coating  
exclusively for ZrO, **shank Ø 6 mm**

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3,0

## Diamant-Instrumente

Diamond coated instruments  
Instruments diamantés  
Strumenti diamantati

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### Zircon oxide abrasives

FG shank, to be used only with spray cooling.

acurata designed the "zircon-a-diamond" FG instruments with a special fine diamond bonding especially tuned for sintered zircon oxide ceramics, to avoid micro cracks and overheating. The "zircon-a-diamond" is proved by  $\text{ZrO}_2$ -expert users.

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**Ultra fine:** 9 µm

Super finishing of composites / compomeres and ceramics.

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**Extra fine:** 25 µm

Finishing of composites, ceramics and filling materials.

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**Fine:** 40–54 µm

Finishing of prepared tooth surfaces and restorations.

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**Medium:** 90–126 µm

Standard applications cavity / crown / in- onlay preparation.

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**Coarse:** 154 µm

Rapid reduction of tooth substance and removal of old fillings.

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**Super coarse:** 180–200 µm

Coarse pregrinding / removal of old fillings.

Separation of none-precious alloys crowns and bridges.

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The total length of the instruments may vary according to length of working part and neck length!

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Flexible disc, double side coated rim

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Flexible disc with saws, double side coated rim

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22,0 / 0,15

## 379

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Disc for plaster with saws, double side coated

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Disc for plaster perforated with saws, double side coated rim

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Disc for plaster, perforated, double side coated rim

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
#### Smart Grinder

Universal grinding tools for all ceramics

Very good edge stability, long service life, Diamond grain in a metal free bond, marginal heat development, smoothly grinding, all-purpose  
Recommended Speed: 10.000–15.000 min<sup>-1</sup>

	855 104 014 524 050	855 104 107 524 050	855 104 173 524 035	855 104 197 524 040	855 104 243 524 035	855 104 372 524 120
Ø mm / L	5,0/5	5,0/12,0	3,5/10,5	4,0/8,0	3,5/6	12,0/2,0

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#### Smart Grinder "MF"

medium grain size, smooth grinding, efficient cutting, smooth surfaces on all ceramic materials – max. speed 12.000 min<sup>-1</sup>. Work with a low contact pressure.

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**Smart Grinder "MC"**

Coarse grain size, fast cutting, fine surfaces on all ceramic materials – maximum speeds 12.000 min<sup>-1</sup>. Work with a low contact pressure.

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**Smart Grinder "SC"**

Very coarse grit size, fast cutting, rapid removal – maximum speeds 12.000 min<sup>-1</sup>. Work with a low contact pressure.

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**Smart Grinder "MAX"**

for all ceramics and hard metals

Easy polishable surfaces, rapidly grinding, coarse grit, ideal for pre-grinding of retention bars, trimming sprues, high edge stability, long service life, particularly economic, diamond grain in a metal free bond, marginal heat development, smooth grinding – Max. speed 12.000 min<sup>-1</sup>. Work with a low contact pressure.

	855 104 014 MAX 180	855 104 304 MAX 220	855 104 371 MAX 220	855 104 372 MAX 220
Ø mm / L	18,0 / 5,0	22,0 / 3,0	22,0 / 1,0	22,0 / 3,0

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**Ceramic polishing, highshine**

Interspersed with diamonds in high concentration. Without plastic center.

	803 104 304 504 260	803 104 304 514 260	803 104 304 524 260	803 104 304 534 260	803 104 371 504 220	803 104 371 514 220
Ø mm / L	26,0 / 2,0	26,0 / 2,0	26,0 / 2,0	26,0 / 2,0	22,0 / 1,0	22,0 / 1,0

	803 104 243 504 030	803 104 243 514 030	803 104 243 524 030	803 104 243 504 045	803 104 243 514 045	803 104 243 524 045
Ø mm / L	3,0 / 7,5	3,0 / 7,5	3,0 / 7,5	4,5 / 13,0	4,5 / 13,0	4,5 / 13,0

	803 104 372 504 110	803 104 372 514 110	803 104 372 524 110	803 104 372 504 170	803 104 372 514 170	803 104 372 524 170
Ø mm / L	11,0 / 2,0	11,0 / 2,0	11,0 / 2,0	17,0 / 3,0	17,0 / 3,0	17,0 / 3,0

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segmented, disc-shaped ceramic polisher in 3 different grits, for all ceramic materials such as ZrO, lithium-silicates or veneering ceramic, flexibly adapts to almost any surface, easy to use, durable, 100% diamond grain. Max. speed 10.000 min<sup>-1</sup>, opt. speed 8.000–10.000 min<sup>-1</sup>. Supple by the segmented shape, structure of the surfaces are preserved.

Vorpölieren – mittlere Körnung  
Pre-polishing – medium grit  
Pré-polissage – grain moyen  
Pre-lucidatura – grana media

Pölfern – feine Körnung  
Polishing – fine grit  
Polissage – grain fin  
Lucidatura – grana fine

Hochglanz – extra feine Körnung  
Highshine – extra fine grit  
Poli brillant – grain extra fin  
Lucidatura a specchio – grana extra fine

Vorpölieren – mittlere Körnung  
Pre-polishing – medium grit  
Pré-polissage – grain moyen  
Pre-lucidatura – grana media

Pölfern – feine Körnung  
Polishing – fine grit  
Polissage – grain fin  
Lucidatura – grana fine

Hochglanz – extra feine Körnung  
Highshine – extra fine grit  
Poli brillant – grain extra fin  
Lucidatura a specchio – grana extra fine

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**Diamond polishers****type "527" (coarse) + "517" (fine)**

Two-step polishing system for polishing of lithium disilicate ceramic, silicate-ceramic, leucite-ceramic. Max. speed 10.000 min<sup>-1</sup>, work with a low contact pressure.

803 104 243 527 045	803 104 303 527 220	803 104 372 527 110	803 104 243 517 045	803 104 303 517 220	803 104 372 517 110
Ø mm / L	4,5 / 13,0	22,0 / 1,0	11,0 / 2,0	4,5 / 13,0	22,0 / 1,0

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Corund, brown – grit size medium for modelling compound and hard metal alloys

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Feinschleifen EM, NEM und MG.

Corund, light brown – Ceramic bond, soft for fine grinding of precious alloys, non precious alloys and modelling compound.

Corund, châtain clair – Alliage céramique, doux  
Meulage fin des métaux précieux, des alliages non précieux et moule.Corund, marrone chiaro – Legame ceramico, tenero  
Molatura fine di metalli preziosi, legami non preziosi e calco.**Empfohlene Drehzahl 25.000 - 35.000 min<sup>-1</sup>**

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Corund, pink – Ceramic bond, hard, for fine grinding  
of precious alloys, non precious alloys and modelling compound.

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Corund blue – Ceramic bond, soft, fine for precious  
and reduced metal alloys. Special abrasive for frameworks.

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Silicon carbide, green – ceramic bond, grit size medium  
for ceramics and precious metal alloys.

**Empfohlene Drehzahl** · Recommended speed  
Vitesses recommandées · Numeri di giri raccomandati **15.000 - 20.000 min<sup>-1</sup>**

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Arkansas, white, for fine grinding of composites, precious alloys

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
22,0 / 3,0

6,0 / 22,0

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Polisher pins for ceramic polishing

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Occlusal surface polisher pins / Adapter

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Silicon carbide, magnesit bond, extra coarse, soft heat-reduced grinding of ceramics and metal alloys.

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acurata instruments are intended for the dental technology sector and must only be used by dental technicians or respective experts who, due their training and experience, are familiar with the use of the products.

### Information for use

Use only technical correct drives with a stable ball bearing! Carefully clamp the instruments as deep as possible and assure that the desired speed is reached prior to usage. Please note the maximum allowed engine speed stated on the packaging label. Work with a contact pressure of maximum 2 N. Avoid levering or canting, otherwise there is risk of fracture. Separate blunt, twisted or damaged instruments immediately and do not use them again. Improper usage increases risk, higher wear and inferior work results! Avoid contact with H<sub>2</sub>O<sub>2</sub>. The respective maximum speed is a theoretical value given on the basis of the geometry of the instruments. The choice of the working speed is determined by the material subsequently to be used the type of handpiece used and the contact pressures and is at user's discretion.

#### Maximum speeds – table of reference values

##### **FG instruments (ISO 314)** ISO Ø:

005–016 max. 450.000 min<sup>-1</sup> • 018–021 max. 300.000 min<sup>-1</sup> • 023–031 max. 160.000 min<sup>-1</sup> • 033–055 max. 120.000 min<sup>-1</sup>

##### **HP and RA instruments (ISO 104 / 204)** ISO Ø:

003–023 max. 50.000 min<sup>-1</sup> • 025–040 max. 40.000 min<sup>-1</sup> • 045–060 max. 30.000 min<sup>-1</sup> • 060–070 max. 25.000 min<sup>-1</sup> • 070–250 max. 25.000 min<sup>-1</sup> • 300–380 max. 20.000 min<sup>-1</sup> • 450 max. 15.000 min<sup>-1</sup>

##### **Polishing instruments** all ISO Ø: max. 5.000–10.000 min<sup>-1</sup>

**Optimum speeds for all instruments depend on the material to be processed. As a general rule 40 – 50% of the maximum instruction is recommended!**

### Important notice for sintered diamond instruments

Clean sintered diamond instruments by using a wet clearing stone No. 807 000 750 before first use. Used and contaminated tools can be restored likewise. During treatment bond abrasion may cause discoloration on ceramics. However the hard FeMn bond and high diamond density will minimize discoloration to a great extent. We recommend a final cleaning of the surface by sandblasting or brief grinding with silicon-carbide trimmers (green, ceramically bonded).

### Important notice for diamond discs

**Minimum speed:** all ISO-diameters min. 10.000 min<sup>-1</sup>

**Maximum speeds:** ISO-diameters 300–380 max. 20.000 min<sup>-1</sup>  
ISO-diameters 450 max. 15.000 min<sup>-1</sup>

Apply little pressure=1N. Wear protective eyewear!!

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## High performance dental instruments to the highest technological level.

Dentists and dental laboratories now value our practically-oriented, needs-oriented products all over the world. We also listen and ask questions – at the same level as our customers. It is in this context that the highly focused, dedicated perfectionists at acurata manufacture absolutely reliable precision instruments which enable ultra-accurate work in the laboratory and in the dental practice. Reliability in the context of a sensitive collaboration between all of our talents – at the technical and human level.

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