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Manufacturer's information for the recycling of sterilizing instruments according to the RKI-directive and DIN EN ISO 17664 - infection prevention in dentistry requirements of hygiene - medical products semi critical B and critical B

Manufacturer

acurata G+K Mahnhardt Dental e. K.
Schulstraße 25, 94169 Thurmansbang
Tel.: +49 (0) 8504 91 17 0
Fax.: +49 (0) 8504 91 17 90
info@acurata-dental.de
www.acurata-dental.de

Products

This manufacturer's information applies to all instruments delivered by acurata G+K Mahnhardt Dental e. K. which are applied for surgical, Periodontic or endodontic measures.

These are polishers, carbide and diamond instruments as well as instruments made of stainless high quality steel as well as root canal instruments Made of steel or nickel titanium. You must prepare before use all unsterile delivered instruments.

Limitation of recycling

The end of the products life cycle is basically determined by the abrasion and damage caused by use. Numerous recycling has no influence on the performance of these instruments.

Storage and transport

Straight after use on the patient put the instruments into a container filled with the applicable cleansing material/ disinfectant. The transport of the instruments to the place of preparation should be made in a special container. For root canal instruments special interim-containers which are fitted with a layer saturated with disinfectant are applicable.

Cleaning and disinfection

According to the directive of the Robert Koch Institute (RKI) it is preferential that further preparation occurs mechanically. Silicone stoppers have to be removed before the preparation of root canal instruments.

Mechanical cleaning

1.1 Applied equipment

- purifier / disinfectant equipment (RDG)
- instrument support for rotary instruments and root canal instruments

1.2 Preparation

Straight after the mechanical preparation of the instruments put them out of the container and clean them under running water. Put the instruments into an appropriate tray. Clean mechanical equipment according to the manufacturer's instructions. Sight check for intactness and cleanness with an appropriate magnification instrument (from experience 8 -10 magnification makes an optic sight check possible). If there are areas of contamination visible after mechanical preparation, repeat cleaning and disinfection until there is no more contamination visible.

2. Alternative manual cleaning

2.1 Used equipment

- cleaning brush
- appropriate detergent or disinfectant for dental instruments with detected disinfectant effect
- ultrasonic bath

2.2 Preparation

Take the instrument out of the container and respectively the interim-tray and clean the surface contamination under running water from the instrument. Remove adhesive contamination completely by constantly turning the instrument under the brush. Clean the instrument under running water. Put the instruments into a strainer element and into the ultrasonic unit filled with detergent or disinfectant. Clean and disinfect according to the directions of the manufacturers of the ultrasonic bath and the detergent and disinfectant. At the end of the residence time clean the instrument with appropriate water (to avoid remnants preferably use fully desalted water (VE)). Dry the instruments preferably with compressed-air. Sight check for intactness and cleanness with an appropriate magnification instrument (from experience 8 - 10 magnification makes an optic sight check possible). If there are areas of contamination visible, repeat cleaning and chemical disinfection until there is no more contamination visible.

2.3 Control and functional check

- Instruments with the following defects are to be separated immediately
- missing diamond bonding (bare spots)
 - blunt and broken out blades
 - shape damages (e.g. twisted instruments, buckled, twisted or broken working part)
 - corroded surfaces

3. Closing steps of preparation

3.1 Medical products semi critical B: thermal disinfection

Final thermal disinfection in the steam sterilizer and store unpacked in appropriate support or strainer tray.

3.2 Medical products critical B: sterile packaging and sterilization

3.2.1 Packaging

There is a choice of the appropriate packaging for the instrument and the sterilization process. Individually packaging: The packaging has to be large enough that the seal is not under tension. Set: Sort the instruments into the provided tray or lay them onto an all-purpose sterilization tray. The instruments must be protected. For packing the tray apply an appropriate process. Instruments with a limited frequency of use should be marked accordingly.

3.2 Sterilisation

Steam sterilisation with the vacuum process at 134° C in a device according to EN 285 to avoid staining and corrosion, the steam has to be free of components. The advised limit of components for drinking water and steam condensate is fixed by EN 285. Consider manufacturer's information.

3.2.3 Transport and storage

The transport and storage of the packaged sterile goods is effected in a dust, humidity and contamination protected manner.

4. Basic annotation

Consider the valid legal regulations of your country concerning the recycling of medical products (e.g. www.rki.de). The manufacturer assures that the processes of recycling explicated above are appropriate for the preparation of the mentioned groups of instruments for their reutilization.

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